Wide range of Nd:YAG Accessories



R35-Nx

Nd:YAG Spot size 2-20 mm MatrixView®



MarcCo S

Nd:YAG Spot size 10 mm Collimated



MarcCo M

Nd:YAG Spot size 24 mm Collimated



MarcCo L

Nd:YAG Spot size 43 mm Collimated



S11 scanner

Nd:YAG Spot size 3,6,9 mm Area 65 mm x 65 mm Frac3®, VERSA



L-runner Pro

Nd:YAG Spot size 3,6,9,11 mm Area 65 mm x 65 mm Frac3®, VERSA, Fotona PIANO® MatrixView® S



Nx-Runner

Nd:YAG Spot sizes 6,9,11 mm Area 65mm x 65 mm Frac3, VERSA, Fotona PIANO® CoolMist™ MatrixView® S



R27

Nd:YAG Surgical HP For 600, 1000 micron fibers



R27C


Nd:YAG Small surgical HP For 600, 1000 micron fibers



For an algorisar or eliable. In a edible and ansistent by good. Ils ethe nathine. It's an absolutely brilliant worknown seandhas nade much siness what ist oday.


R35-Nx Handpiece

 Integrated CoolMist[™] cooling technology* and MatrixView® thermal sensor




* For a greater cooling effect, a special cold air cooling adapter is available. It can also be used without dry spray.

Nd:YAG Scanner



Er:YAG Accessories



Dynamic Vacuum Technologie

- The system's integrated vacuum technology is based on dynamically controlled intermittent suction and represents a unique add-on to enhance Nx line's exceptional versatility.
- Works with a variety of radial applicators for highly efficient deep-tissue body and facial aesthetic treatments.

Body S Body M Face S Face XS

