

STARFORMER

Magnetic Body Sculpting and Muscle Strengthening



Fotona[®]
choose perfection

STARFORMER




StarFormer is Fotona's latest innovation in body shaping. This breakthrough technology utilizes magnetic pulse energy to stimulate muscle contraction. As a result, muscles strengthen, tone and firm.

The system provides two treatment options:

- **TightWave™**
muscle strengthening, toning and firming of weakened areas
- **IntimaWave™**
strengthening of pelvic floor muscles to treat stress urinary incontinence

How it Works

StarFormer uses magnetic stimulation technology to selectively target muscle tissues. Magnetic pulses harmlessly penetrate into the muscle layers causing them to automatically contract. The end result is increased muscle mass.



TIGHTWAVE™

TightWave™ is a breakthrough treatment for effective strengthening, shaping and firming of 7 different areas of the body:



- Glutes
- Abdomen
- Hamstring & Thighs
- Calves
- Biceps & Triceps
- Shoulder
- Back




The treatment is performed with the patient lying down fully clothed. One session generally lasts between 20 to 30 minutes with patients resuming their daily activities immediately after treatment.

INTIMAWAVE™

StarFormer can additionally be used to stimulate pelvic floor muscles and back, utilizing two applicators to treat all types of male or female urinary and faecal incontinence. During treatment the patient remains fully clothed and sits on a comfortable IntimaWave™ chair. The magnetic applicators in the chair stimulate muscle contraction in the entire pelvic floor area, increasing the strength and endurance of the pelvic floor. The result is effective and long-term relief of uncontrolled urine flow.




Before




Pelvic floor muscles insufficiently support pelvic organs effecting bladder control.

After




Magnetic stimulation induces pelvic floor muscle strengthening and the return of normal continence function.



One Device, Multiple Solutions

Key Benefits

- Safe, quick and easy muscle strengthening
- Non-invasive
- Improvement after only a few sessions
- Walk-in/walk-out procedure
- High success rate and patient satisfaction



| | |
|-------------------------------|--------------------------|
| Device name: | Fotona StarFormer |
| Input voltage: | 100-240V AC, 50~60Hz |
| Power: | 1500 VA |
| IP protection: | IPX0 |
| Size (W x L x H): | 421 mm x 630 mm x 843 mm |
| Weight: | 50 kg |
| IEC 60601/1 classification: | Class I, Type BF |
| MDD 93/42/EEC classification: | Class IIa |

www.fotona.com

Fotona, d. o. o.
Stegne 7
1000 Ljubljana
Slovenia
EU

Fotona, LLC
2307 Springlake Road
#518
Dallas, TX 75234
USA

Fotona GmbH
Hohlbachweg 2
73344 Gruibingen,
Germany
EU

Fotona Beauty Light, (Suzhou)
Medical Devices Co, Ltd.
No 2, Zengfu Road, Guli Town
Changshu City, Jiangsu Province
CHINA, 215515



2265 For related patents see: www.fotona.com/patents